

Invacare® | **Birdie® Family**

NEW



Birdie



Birdie Compact



Birdie Plus

Transportability and patient comfort go hand-in-hand

With maximum space for the user, the *Invacare Birdie*, *Birdie Compact* and *Birdie Plus* offer a truly comfortable lift and transfer to or from beds, chairs or even the floor. The lifters are designed to ensure that folding and unfolding can be carried out easily and without the need for tools. In addition, the lifter can be dismantled into two parts if needed – also without tools.



Yes, you can.®



Total user comfort

Maximizes space for the user even when the hoist is in the highest position. This gives a very comfortable lift and transfer since the user can easily be rotated in the hoist. Consequently, it makes the placing of the user into the bed/wheelchair much easier for the carer. The optimization of space creates a high level of comfort since the user's knees will not be close to the actuator of the hoist. You have, then, a very low risk of injuring the knees of the user.

Ease of handling: simple to fold, unfold and dismantle without tools

Handling the lifter has been made as smooth and easy as possible as the hoists can easily be folded and unfolded without the need for tools. Hoists take up minimal space when stored and are easy to push and transport. If required, the mast and boom can be easily detached from the base (no tools required).

TÜV-certification for maximum security

For total security, the hoists come with a snap-lock system, thus minimizing the risk of accident. The lifters are CE-certified in line with the Medical Device Directive 93/42/ECC. To ensure

the hoists comply with EN ISO 10535 (the European standard for hoists) the hoist has been rigorously tested both by Invacare and the German TÜV institute, an accredited external body.

Features and Options



Long boom – high user comfort

The user can easily be rotated and positioned. No risk of injuries to the knees.



Easy to fold and dismantle

Ease of transport and storage. No tools required when folding / unfolding or when dismantling the hoist.



Snap lock

Optimal security to avoid risk of accidents.



Emergency stop, emergency lowering and integrated charger

Manual or electric emergency lowering (compact version features only electric lowering). To charge the hoist, simply plug in.

Features and Options



Only available on Birdie Plus:

- Even Easier manoeuvring - 125mm rear wheels
- Electrical emergency lifting
- Showing battery status
- Soft start and stop
- Covered electrical legspread - noise reduction



Click on spreader bar

Easy to change the spreader bar.



Manual leg spread

Simply step on the pedal to adjust the leg spread.



Handle for leg spread

Handle can easily be attached to operate the manual leg spread.



Electrical leg spread

Available for Invacare **Birdie** Standard version.



2 point quick release spreader bars

Easy to attach and detach. Available widths: 35 cm, 45 cm, 55 cm.



Push bar and hand control

Superior ergonomics for carers.



Castors

100 mm or 75 mm castors available (braking castors are shown). **Birdie Plus** comes with 100 mm front castors and 125 mm rear castors.



Scale






If needed the scale is available approved, according to the medical standard EN 45501, certified and sealed in order to be used to medicate.








4 point quick release spreader bars

Easy to attach and detach. Available widths: 45 cm and 55 cm.

Technical data

					
Birdie	115 cm	65,5 cm	42 kg	46,5 - 168,5 cm	10 cm
Birdie Compact	100 cm	53 cm	36 kg	51,5 - 158,5 cm	10 cm
Birdie Plus	116,5 cm	68 cm	44 kg	47,5 - 170 cm	11,5 cm

					
Birdie	104 cm	170 kg	40 cycles	2,9 Ah	140 cm
Birdie Compact	89 cm	150 kg	40 cycles	2,9 Ah	107 cm
Birdie Plus	104 cm	170 kg	40 cycles	2,9 Ah	140 cm

All measures are taken with 75 mm castors and 2 point spreaderbar.
For 100 mm castors add 1,5 cm for height measures and 2 cm for width.

Frame colour



Light grey /
Dark grey

Invacare International Sarl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.

Birdie Family - EU - 08/2010

